



DONJOY®




T 800.336.6569 D 760.727.1280 F 800.936.6569
2900 Lake Vista Drive | Dallas, TX 75067 | U.S.A.
DJOGlobal.com/donjoy


Copyright © 2019 by DJO, LLC
MKT00-8551 Rev A

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

iceman®
CLASSIC³



ICEMAN with dual pump recirculation technology. The DonJoy® IceMan® features a patented recirculation system that allows water warmed after flowing through the cold pad to be preserved and remixed with cooler ice water at a constant flow rate, providing consistent cool water distribution throughout the cold pad.



COMPETITORS Other versions restrict or pinch water flow to the cold pad, causing variable flow rates and inconsistent cold pad temperatures.



WARNING!

THIS DEVICE CAN BE COLD ENOUGH TO CAUSE SERIOUS INJURY. Please read and understand all warnings and instructions for use before using the device. DO NOT use this device without a prescription and direction from a physician.

DONJOY® ICEMAN CLASSIC3™

The IceMan CLASSIC3 cold therapy unit helps reduce pain and swelling, speeding up rehabilitation and recovery. The IceMan® helps provide extended cold therapy to aid a variety of indications and protocols as directed by a medical professional. It utilizes DonJoy's patented recirculation system, which helps maintain more consistent and accurate temperatures than other competitive cold therapy units, in a pre-set configuration.

PROVEN PERFORMANCE

- Time tested with over 20-years of experience delivering industry leading cold therapy solutions
- Developed through direct customer interaction and feedback

TECHNOLOGY ADVANTAGE

- Dual pump recirculation system helps deliver consistent cold therapy throughout the cold pad, eliminating freezing inlet temperatures

PATIENT CONVENIENCE AND EASE-OF-USE

- Self-priming, plug-n-play operation
- Designed for easy portability and storage
- A complete line of cold pads to deliver cold therapy to where it is needed the most

The use of cold therapy may reduce the need for narcotics and help accelerate rehabilitation and recovery.

STEPS:

1. Unlock lid by rotating the handle forward.
2. Add ice to fill line inside the device.
3. Add cold water to fill line.
4. With handle in forward position place lid on device. Press lid down while rotating handle up and back to latch closed.
5. Connect the IceMan CLASSIC3 hose to the cold pad hose. To ensure a reliable connection, "snap" or "click" hoses together into place so that the fit is tight and snug.
6. To turn device on, insert cord into the connection port on the back of the device and plug the power supply into the wall outlet.